

兔抗 PDP1 多克隆抗体

中文名称：兔抗 PDP1 多克隆抗体

英文名称：Anti-PDP1 rabbit polyclonal antibody

别名：PDH; PDP; PDPC; PPM2C

相关类别：一抗

储存：冷冻（-20℃）

宿主：Rabbit

抗原：PDP1

反应种属：Human, Mouse, Rat

标记物：Unconjugate

克隆类型：rabbit polyclonal

技术规格

Background:	Pyruvate dehydrogenase (E1) is one of the three components (E1, E2, and E3) of the large pyruvate dehydrogenase complex. Pyruvate dehydrogenase kinases catalyze phosphorylation of serine residues of E1 to inactivate the E1 component and inhibit the complex. Pyruvate dehydrogenase phosphatases catalyze the dephosphorylation and activation of the E1 component to reverse the effects of pyruvate dehydrogenase kinases.
Applications:	ELISA, IHC
Name of antibody:	PDP1
Immunogen:	Fusion protein of human PDP1

Full name:	pyruvate dehydrogenase phosphatase catalytic subunit 1
Synonyms :	PDH; PDP; PDPC; PPM2C
SwissProt:	Q9P0J1
ELISA Recommended dilution:	2000-5000
IHC positive control:	Human thyroid cancer and human cervical cancer
IHC Recommend dilution:	30-150

