

兔抗 PER2 多克隆抗体

中文名称：兔抗 PER2 多克隆抗体

英文名称：Anti-PER2 rabbit polyclonal antibody

别名：FASPS; FASPS1

相关类别：一抗

储存：冷冻（-20℃）

宿主：Rabbit

抗原：PER2

反应种属：Human

标记物：Unconjugate

克隆类型：rabbit polyclonal

技术规格

Background:	This gene is a member of the Period family of genes and is expressed in a circadian pattern in the suprachiasmatic Nucleus, the primary circadian pacemaker in the mammalian brain. Genes in this family encode components of the circadian rhythms of locomotor activity, metabolism, and behavior. This gene is upregulated by CLOCK/ARNTL heterodimers but then represses this upregulation in a feedback loop using PER/CRY heterodimers to interact with CLOCK/ARNTL. Polymorphisms in this gene may increase the risk of getting certain cancers and have been linked to sleep disorders.
Applications:	ELISA, IHC
Name of antibody:	PER2

Immunogen:	Synthetic peptide of human PER2
Full name:	period circadian clock 2
Synonyms :	FASPS; FASPS1
SwissProt:	O15055
ELISA Recommended dilution:	2000-10000
IHC positive control:	Human thyroid cancer and Human lung cancer
IHC Recommend dilution:	100-300

