

兔抗 WDSUB1 多克隆抗体

中文名称：兔抗 WDSUB1 多克隆抗体

英文名称：Anti-WDSUB1 rabbit polyclonal antibody

别名：UBOX6; WDSAM1

相关类别：一抗

储存：冷冻（-20℃）

宿主：Rabbit

抗原：WDSUB1

反应种属：Human

标记物：Unconjugate

克隆类型：rabbit polyclonal

技术规格

Background:

WDSUB1 (WD repeat, SAM and U-box domain-containing protein 1), also known as UBOX6 or WDSAM1, is a 476 amino acid protein that contains one SAM (sterile alpha motif) domain, one U-box domain and seven WD repeats. Existing as two isoforms due to alternative splicing, WDSUB1 is encoded by a gene located on chromosome 2. The second largest human chromosome, chromosome 2 encodes over 1,400 genes and comprises nearly 8% of the human genome, housing a number of disease-associated genes. Harlequin ichthyosis, a rare and morbid skin deformity, is associa

	ted with mutations in the ABCA12 gene, while the lipid metabolic disorder sitosterolemia is associated with defects in the ABCG5 and ABCG8 genes. Additionally, an extremely rare recessive genetic disorder, Alström syndrome, is caused by mutations in the ALMS1 gene, which maps to chromosome 2.
Applications:	ELISA, WB, IHC
Name of antibody:	WDSUB1
Immunogen:	Synthetic peptide of human WDSUB1
Full name:	WD repeat, sterile alpha motif and U-box domain containing 1
Synonyms:	UBOX6; WDSAM1
SwissProt:	Q8N9V3
ELISA Recommended dilution:	5000-10000
IHC positive control:	Human liver cancer
IHC Recommend dilution:	25-100
WB Predicted band size:	53 kDa
WB Positive control:	231 , K562 , 293T and Hela cell lysates
WB Recommended dilution:	200-1000



