

兔抗 KBTBD8 多克隆抗体

- 中文名称：兔抗 KBTBD8 多克隆抗体
- 英文名称：Anti-KBTBD8 rabbit polyclonal antibody
- 别名：kelch repeat and BTB domain containing 8; TAKRP; TA-KRP
- 相关类别：一抗
- 储存：冷冻（-20℃）
- 宿主：Rabbit
- 抗原：KBTBD8
- 反应种属：Human, Mouse, Rat
- 标记物：Unconjugate
- 克隆类型：rabbit polyclonal

技术规格

Background:	Substrate-specific adapter of a BCR (BTB-CUL3-RBX1) E3 ubiquitin ligase complex that acts as a regulator of neural crest specification. The BCR(KBTBD8) complex acts by mediating monoubiquitination of NOLC1 and TCOF1: monoubiquitination promotes the formation of a NOLC1-TCOF1 complex that acts as a platform to connect RNA polymerase I with enzymes responsible for ribosomal processing and modification, leading to remodel the translational program of differentiating cells in favor of neural crest specification.
Applications:	ELISA, WB
Name of antibody:	KBTBD8

Immunogen:	Fusion protein of human KBTBD8
Full name:	kelch repeat and BTB domain containing 8
Synonyms:	TAKRP; TA-KRP
SwissProt:	Q8NFI9
ELISA Recommended dilution:	5000-10000
WB Predicted band size:	69 kDa
WB Positive control:	NIH/3T3 cell lysate
WB Recommended dilution:	1000-5000

