

兔抗 KBTBD11 多克隆抗体

中文名称：兔抗 KBTBD11 多克隆抗体

英文名称：Anti-KBTBD11 rabbit polyclonal antibody

别名：kelch repeat and BTB domain containing 11; KLHDC7C

相关类别：一抗

储存：冷冻（-20℃）

宿主：Rabbit

抗原：KBTBD11

反应种属：Human, Mouse

标记物：Unconjugate

克隆类型：rabbit polyclonal

技术规格

Background:	The function of this protein remains unknown.
Applications:	ELISA, WB
Name of antibody:	KBTBD11
Immunogen:	Synthetic peptide of human KBTBD11
Full name:	kelch repeat and BTB domain containing 11
Synonyms:	KLHDC7C
SwissProt:	O94819
ELISA Recommended dilution:	5000-10000
WB Predicted band size:	66 kDa
WB Positive control:	Human metastatic papillary carcinoma(thyroid cancer) tissue I

	ysate
WB Recommended dilution :	1000-5000

