

## 兔抗 SGK1 多克隆抗体

中文名称：兔抗 SGK1 多克隆抗体

英文名称：Anti-SGK1 rabbit polyclonal antibody

别名：SGK

相关类别：一抗

储存：冷冻（-20℃）

宿主：Rabbit

抗原：SGK1

反应种属：Human, Mouse, Rat

标记物：Unconjugate

克隆类型：rabbit polyclonal

### 技术规格

<b>Background:</b>	This gene encodes a serine/threonine protein kinase that plays an important role in cellular stress response. This kinase activates certain potassium, sodium, and chloride channels, suggesting an involvement in the regulation of processes such as cell survival, neuronal excitability, and renal sodium excretion. High levels of expression of this gene may contribute to conditions such as hypertension and diabetic nephropathy. Several alternatively spliced transcript variants encoding different isoforms have been noted for this gene.
<b>Applications:</b>	ELISA, IHC

<b>Name of antibody:</b>	SGK1
<b>Immunogen:</b>	Synthetic peptide of human SGK1
<b>Full name:</b>	serum/glucocorticoid regulated kinase 1
<b>Synonyms :</b>	SGK
<b>SwissProt:</b>	O00141
<b>ELISA Recommended dilution:</b>	1000-2000
<b>IHC positive control:</b>	Human liver cancer and Human lung cancer
<b>IHC Recommend dilution:</b>	25-100

