

兔抗 SH3KBP1 多克隆抗体

中文名称: 兔抗 SH3KBP1 多克隆抗体

英文名称: Anti-SH3KBP1 rabbit ployclonal antibody

别 名: HSB1; CIN85; GIG10; HSB-1; MIG18; CD2BP3

相关类别: 一抗

储 存: 冷冻(-20℃)

宿 主: Rabbit

抗原: SH3KBP1

反应种属: Human, Mouse, Rat

标记物: Unconjugate

克隆类型: rabbit polyclonal

技术规格

Background:	This gene encodes an adapter protein that contains three N-terminal Src homology domains, a proline rich region a nd a C-terminal coiled-coil domain. The encoded protein facilitates protein-protein interactions and has been implic ated in numerous cellular processes including apoptosis, c ytoskeletal rearrangement, cell adhesion and in the regula tion of clathrin-dependent endocytosis. Alternate splicing results in multiple transcript variants.
Applications:	ELISA, WB, IHC
Name of antibody:	SH3KBP1
Immunogen:	Fusion protein of human SH3KBP1



Full name:	SH3-domain kinase binding protein 1
Synonyms:	HSB1; CIN85; GIG10; HSB-1; MIG18; CD2BP3
SwissProt:	Q96B97
ELISA Recommended dilution:	2000-5000
IHC positive control:	Human liver cancer and human breast cancer
IHC Recommend dilution:	30-150
WB Predicted band size:	73 kDa
WB Positive control:	Jurkat cells
WB Recommended dilution:	500-2000





