

兔抗 SHBG 多克隆抗体

中文名称：兔抗 SHBG 多克隆抗体

英文名称： Anti-SHBG rabbit polyclonal antibody

别名： ABP, SBP, TEBG

相关类别： 一抗

储存： 冷冻（-20℃）

宿主： Rabbit

抗原： SHBG

反应种属： Human, Mouse, Rat

标记物： Unconjugate

克隆类型： rabbit polyclonal

技术规格

Background:	This gene encodes a steroid binding protein that was first described as a plasma protein secreted by the liver but is now thought to participate in the regulation of steroid responses. The encoded protein binds each steroid molecule as a dimer formed from identical or nearly identical monomers. The use of alternate promoters and alternatively spliced transcripts have been described. Multiple transcript variants encoding different isoforms have been found for this gene.
Applications:	ELISA, IHC

Name of antibody:	SHBG
Immunogen:	Synthetic peptide of human SHBG
Full name:	Sex hormone-binding globulin
Synonyms :	ABP, SBP, TEBG
SwissProt:	P04278
ELISA Recommended dilution:	1000-5000
IHC positive control:	Human colorectal cancer and Human ovarian cancer
IHC Recommend dilution:	25-100

