

兔抗 PIWIL3 多克隆抗体

中文名称：兔抗 PIWIL3 多克隆抗体

英文名称：Anti-PIWIL3 rabbit polyclonal antibody

别名：HIWIL3

相关类别：一抗

储存：冷冻（-20℃）

抗原：PIWIL3

宿主：Rabbit

反应种属：Human

标记物：Unconjugate

克隆类型：rabbit polyclonal

技术规格

Background:	This gene encodes a member of the PIWI subfamily of Argonaute family proteins. This subfamily of proteins contains a PAZ domain, found in proteins involved in RNA-mediated gene silencing, and a C-terminal Piwi domain. The encoded protein is thought to function in maintenance of germline cells. Multiple transcript variants encoding different isoforms have been found for this gene.
Applications:	ELISA, IHC
Name of antibody:	PIWIL3
Immunogen:	Synthetic peptide of human PIWIL3
Full name:	piwi like RNA-mediated gene silencing 3

Synonyms:	HIWI3
SwissProt:	Q7Z3Z3
ELISA Recommended dilution:	5000-10000
IHC positive control:	Human liver cancer and Human ovarian cancer
IHC Recommend dilution:	50-100



