

兔抗 PLIN2 多克隆抗体

中文名称：兔抗 PLIN2 多克隆抗体

英文名称：Anti-PLIN2 rabbit polyclonal antibody

别名：ADFP; ADRP

相关类别：一抗

储存：冷冻（-20℃）

宿主：Rabbit

抗原：PLIN2

反应种属：Human, Mouse

标记物：Unconjugate

克隆类型：rabbit polyclonal

技术规格

Background:	The protein encoded by this gene belongs to the perilipin family, members of which coat intracellular lipid storage droplets. This protein is associated with the lipid globule surface membrane material, and maybe involved in development and maintenance of adipose tissue. However, it is not restricted to adipocytes as previously thought, but is found in a wide range of cultured cell lines, including fibroblasts, endothelial and epithelial cells, and tissues, such as lactating mammary gland, adrenal cortex, Sertoli and Leydig cells, and hepatocytes in alcoholic liver cirrhosis, suggesting that it may serve as a marker of lipid accumu
--------------------	---

	lation in diverse cell types and diseases.
Applications:	ELISA, WB, IHC
Name of antibody:	PLIN2
Immunogen:	Synthetic peptide of human PLIN2
Full name:	perilipin 2
Synonyms :	ADFP; ADRP
SwissProt:	Q99541
ELISA Recommended dilution:	2000-5000
IHC positive control:	Human thyroid cancer and Human tonsil
IHC Recommend dilution:	50-200
WB Predicted band size:	48 kDa
WB Positive control:	Mouse liver tissue
WB Recommended dilution:	500-2000



