

兔抗 KIAA0513 多克隆抗体

中文名称：兔抗 KIAA0513 多克隆抗体

英文名称：Anti-KIAA0513 rabbit polyclonal antibody

别名：KIAA0513

相关类别：一抗

储存：冷冻（-20℃）

宿主：Rabbit

抗原：KIAA0513

反应种属：Human, Mouse

标记物：Unconjugate

克隆类型：rabbit polyclonal

技术规格

Background:	The function of this protein remains unknown.
Applications:	ELISA, IHC
Name of antibody:	KIAA0513
Immunogen:	Fusion protein of human KIAA0513
Full name:	KIAA0513
SwissProt:	O60268
ELISA Recommended dilution:	5000-10000
IHC positive control:	Human tonsil and Human brain
IHC Recommend dilution:	50-300

