

兔抗 SLC7A4 多克隆抗体

中文名称：兔抗 SLC7A4 多克隆抗体

英文名称：Anti-SLC7A4 rabbit polyclonal antibody

别名：VH; CAT4; CAT-4; HCAT3

抗原：SLC7A4

储存：冷冻（-20℃）避光

宿主：Rabbit

反应种属：Human

相关类别：一抗

标记物：Unconjugate

克隆类型：Unconjugate

技术规格

Background:	Involved in the transport of the cationic amino acids (arginine, lysine and ornithine).
Applications:	WB
Name of antibody:	SLC7A4
Immunogen:	Synthesized peptide derived from C-terminal of human SLC7A4.
Full name:	solute carrier family 7, member 4
Synonyms :	VH; CAT4; CAT-4; HCAT3
SwissProt:	O43246
WB Predicted band size:	68 kDa
WB Positive control:	LOVO cells and A549 cells lysates

WB Recommended dilution: 500-3000

