

兔抗 SLC12A6 多克隆抗体

中文名称：兔抗 SLC12A6 多克隆抗体

英文名称：Anti-SLC12A6 rabbit polyclonal antibody

别名：KCC3, ACCPN, KCC3A, KCC3B

相关类别：一抗

储存：冷冻（-20℃）

宿主：Rabbit

抗原：SLC12A6

反应种属：Human

标记物：Unconjugate

克隆类型：rabbit polyclonal

技术规格

Background:	This gene is a member of the K-Cl cotransporter (KCC) family. K-Cl cotransporters are integral membrane proteins that lower intracellular chloride concentrations below the electrochemical equilibrium potential. The proteins encoded by this gene are activated by cell swelling induced by hypotonic conditions. Alternate splicing results in multiple transcript variants encoding different isoforms. Mutations in this gene are associated with agenesis of the corpus callosum with peripheral neuropathy.
Applications:	ELISA, IHC

Name of antibody:	SLC12A6
Immunogen:	Synthetic peptide of human SLC12A6
Full name:	Solute carrier family 12 (potassium/chloride transporter s), member 6
Synonyms :	KCC3, ACCPN, KCC3A, KCC3B
SwissProt:	Q9UHW9
ELISA Recommended dilution:	1000-5000
IHC positive control:	Human breast cancer
IHC Recommend dilution:	25-100

