

兔抗 SLC16A12 多克隆抗体

中文名称：兔抗 SLC16A12 多克隆抗体

英文名称： Anti-SLC16A12 rabbit polyclonal antibody

别 名： CJMG; MCT12

相关类别： 一抗

储 存： 冷冻（-20℃）

宿 主： Rabbit

抗 原： SLC16A12

反应种属： Human

标记物： Unconjugate

克隆类型： rabbit polyclonal

技术规格

Background:	This gene encodes a transmembrane transporter that likely plays a role in monocarboxylic acid transport. A mutation in this gene has been associated with juvenile cataracts with microcornea and renal glucosuria.
Applications:	ELISA, IHC
Name of antibody:	SLC16A12
Immunogen:	Synthetic peptide of human SLC16A12
Full name:	solute carrier family 16, member 12
Synonyms :	CJMG; MCT12
SwissProt:	Q6ZSM3
ELISA Recommended dilution:	1000-2000

IHC positive control:	Human liver cancer
IHC Recommend dilution:	25-100

