

## 兔抗 SLC16A14 多克隆抗体

中文名称：兔抗 SLC16A14 多克隆抗体

英文名称：Anti-SLC16A14 rabbit polyclonal antibody

别名：MCT14

相关类别：一抗

储存：冷冻（-20℃）

宿主：Rabbit

抗原：SLC16A14

反应种属：Human, Mouse

标记物：Unconjugate

克隆类型：rabbit polyclonal

### 技术规格

|                                    |                                                                                                                       |
|------------------------------------|-----------------------------------------------------------------------------------------------------------------------|
| <b>Background:</b>                 | Proton-linked monocarboxylate transporter. May catalyze the transport of monocarboxylates across the plasma membrane. |
| <b>Applications:</b>               | ELISA, WB, IHC                                                                                                        |
| <b>Name of antibody:</b>           | SLC16A14                                                                                                              |
| <b>Immunogen:</b>                  | Synthetic peptide of human SLC16A14                                                                                   |
| <b>Full name:</b>                  | solute carrier family 16 member 14                                                                                    |
| <b>Synonyms:</b>                   | MCT14                                                                                                                 |
| <b>SwissProt:</b>                  | Q7RTX9                                                                                                                |
| <b>ELISA Recommended dilution:</b> | 5000-10000                                                                                                            |
| <b>IHC positive control:</b>       | Human esophagus cancer and Human colorectal cancer                                                                    |

|                                 |                            |
|---------------------------------|----------------------------|
|                                 | er                         |
| <b>IHC Recommend dilution:</b>  | 25-100                     |
| <b>WB Predicted band size:</b>  | 57 kDa                     |
| <b>WB Positive control:</b>     | Mouse kidney tissue lysate |
| <b>WB Recommended dilution:</b> | 200-1000                   |



