

## 兔抗 SLC17A7 多克隆抗体

中文名称: 兔抗 SLC17A7 多克隆抗体

英文名称: Anti-SLC17A7 rabbit polyclonal antibody

别 名: BNPI; VGLUT1

相关类别: 一抗

储 存: 冷冻(-20℃)

宿 主: Rabbit

抗 原: SLC17A7

反应种属: Human, Mouse, Rat

标 记 物: Unconjugate

克隆类型: rabbit polyclonal

## 技术规格

Background:	The protein encoded by this gene is a vesicle-bound, sodium-dependent phosphate transporter that is specifically expressed in the neuron-rich regions of the brain. It is preferentially associated with the membranes of synaptic vesicles and functions in glutamate transport. The protein shares 82% identity with the differentiation-associated Na-dependent inorganic phosphate cotransporter and they appear to form a distinct class within the Na+/Pi cotransporter family.
Applications:	ELISA, IHC
Name of antibody:	SLC17A7



Immunogen:	Synthetic peptide of human SLC17A7
Full name:	solute carrier family 17 (sodium-dependent inorganic p hosphate cotransporter), member 7
Synonyms:	BNPI; VGLUT1
SwissProt:	Q9P2U7
ELISA Recommended dilution:	5000-240000
IHC positive control:	Human brain
IHC Recommend dilution:	200-500

