

兔抗 KLC2 多克隆抗体

中文名称： 兔抗 KLC2 多克隆抗体

英文名称： Anti-KLC2 rabbit polyclonal antibody

别 名： kinesin light chain 2

相关类别： 一抗

储 存： 冷冻（-20℃）

宿 主： Rabbit

抗 原： KLC2

反应种属： Human, Mouse

标记物： Unconjugate

克隆类型： rabbit polyclonal

技术规格

Background:	The protein encoded by this gene is a light chain of kinesin, a molecular motor responsible for moving vesicles and organelles along microtubules. Defects in this gene are a cause of spastic paraplegia, optic atrophy, and neuropathy (SPOAN) syndrome.
Applications:	ELISA, WB, IHC
Name of antibody:	KLC2
Immunogen:	Fusion protein of human KLC2
Full name:	kinesin light chain 2
SwissProt:	Q9H0B6
ELISA Recommended dilution:	5000-10000

IHC positive control:	Human thyroid cancer
IHC Recommend dilution:	150-300
WB Predicted band size:	69 kDa
WB Positive control:	Human cerebrum tissue and Mouse brain tissue lysates
WB Recommended dilution:	500-2000



