

兔抗 ZIC1 多克隆抗体

中文名称：兔抗 ZIC1 多克隆抗体

英文名称：Anti-ZIC1 rabbit polyclonal antibody

级 别：ZIC; ZNF201

抗 原：ZIC1

储 存：冷冻（-20℃）

宿 主：Rabbit

相关类别：一抗

反应种属：Human, Mouse

标 记 物：Unconjugate

克隆类型：rabbit polyclonal

技术规格

Background:	This gene encodes a member of the ZIC family of C2H2-type zinc finger proteins. Members of this family are important during development. Aberrant expression of this gene is seen in medulloblastoma, a childhood brain tumor. This gene is closely linked to the gene encoding zinc finger protein of the cerebellum 4, a related family member on chromosome 3. This gene encodes a transcription factor that can bind and transactivate the apolipoprotein E gene.
Applications:	ELISA, IHC
Name of antibody:	ZIC1

Immunogen:	Synthetic peptide of human ZIC1
Full name:	Zic family member 1
Synonyms:	ZIC; ZNF201
SwissProt:	Q15915
ELISA Recommended dilution:	5000-10000
IHC positive control:	Human liver cancer and human esophagus cancer
IHC Recommend dilution:	25-100



