

兔抗 SLC30A6 多克隆抗体

- 中文名称：兔抗 SLC30A6 多克隆抗体
- 英文名称：Anti-SLC30A6 rabbit polyclonal antibody
- 别名：solute carrier family 30 member 6; ZNT6; MST103; MSTP103
- 相关类别：一抗
- 储存：冷冻（-20℃）
- 宿主：Rabbit
- 抗原：SLC30A6
- 反应种属：Human, Mouse
- 标记物：Unconjugate
- 克隆类型：rabbit polyclonal

技术规格

Background:	This gene encodes a member of a family of proteins that function as zinc transporters. This protein can regulate subcellular levels of zinc in the Golgi and vesicles. Expression of this gene is altered in the Alzheimer's disease brain plaques.
Applications:	ELISA, IHC
Name of antibody:	SLC30A6
Immunogen:	Fusion protein of human SLC30A6
Full name:	solute carrier family 30 member 6
Synonyms:	ZNT6; MST103; MSTP103
SwissProt:	Q6NXT4

ELISA Recommended dilution:	5000-10000
IHC positive control:	Human thyroid cancer
IHC Recommend dilution:	50-200

