

兔抗 KLRF1 多克隆抗体

中文名称: 兔抗 KLRF1 多克隆抗体

英文名称: Anti-KLRF1 rabbit polyclonal antibody

别名: NKp80; CLEC5C

相关类别: 一抗

储存: 冷冻 (-20℃)

宿主: Rabbit

抗原: KLRF1

反应种属: Human

标记物: Unconjugate

克隆类型: rabbit polyclonal

技术规格

Background:	KLRF1, an activating homodimeric C-type lectin-like receptor (CTLR), is expressed on nearly all natural killer (NK) cells and stimulates their cytotoxicity and cytokine release
Applications:	ELISA, IHC
Name of antibody:	KLRF1
Immunogen:	Synthetic peptide of human KLRF1
Full name:	killer cell lectin-like receptor subfamily F, member 1
Synonyms :	NKp80; CLEC5C
SwissProt:	Q9NZS2

ELISA Recommended dilution:	2000-5000
IHC positive control:	Human ovarian cancer and Human thyroid cancer
IHC Recommend dilution:	50-200

