

兔抗 KRT7 多克隆抗体

中文名称: 兔抗 KRT7 多克隆抗体

英文名称: Anti-KRT7 rabbit polyclonal antibody

别名: K7; CK7; SCL; K2C7

抗原: KRT7

储存: 冷冻 (-20℃)

宿主: Rabbit

反应种属: Human

相关类别: 一抗

标记物: Unconjugate

克隆类型: rabbit polyclonal

技术规格

Background:

The protein encoded by this gene is a member of the keratin gene family. The type II cytokeratins consist of basic or neutral proteins which are arranged in pairs of heterotypic keratin chains coexpressed during differentiation of simple and stratified epithelial tissues. This type II cytokeratin is specifically expressed in the simple epithelia lining the cavities of the internal organs and in the gland ducts and blood vessels. The genes encoding the type II cytokeratins are clustered in a region of chromosome 12q12-q13. Alternative splicing may result in several transcript variants; ho

	wever, not all variants have been fully described.
Applications:	ELISA, WB, IHC
Name of antibody:	KRT7
Immunogen:	Synthetic peptide of human KRT7
Full name:	keratin 7, type II
Synonyms:	K7; CK7; SCL; K2C7
SwissProt:	P08729
IHC positive control:	Human esophagus cancer and human thyroid cancer
IHC Recommend dilution:	10-50
WB Predicted band size:	51 kDa
WB Positive control:	Hela cell lysate
WB Recommended dilution:	200-1000





