

兔抗 ZNF35 多克隆抗体

中文名称：兔抗 ZNF35 多克隆抗体

英文名称：Anti-ZNF35 rabbit polyclonal antibody

别名：zinc finger protein 35; HF10; HF.10; Zfp105

相关类别：一抗

储存：冷冻（-20℃）

宿主：Rabbit

抗原：ZNF35

反应种属：Human

标记物：Unconjugate

克隆类型：rabbit polyclonal

技术规格

Background:	May be involved in transcriptional regulation. Involved in cell differentiation and/or proliferation.
Applications:	ELISA, IHC
Name of antibody:	ZNF35
Immunogen:	Synthetic peptide of human ZNF35
Full name:	zinc finger protein 35
Synonyms:	HF10; HF.10; Zfp105
SwissProt:	P13682
ELISA Recommended dilution:	5000-10000
IHC positive control:	Human liver cancer
IHC Recommend dilution:	50-100

