

兔抗 ZNF43 多克隆抗体

中文名称：兔抗 ZNF43 多克隆抗体

英文名称：Anti-ZNF43 rabbit polyclonal antibody

别名：HTF6; KOX27; ZNF39L1

抗原：ZNF43

相关类别：一抗

储存：冷冻（-20℃）

宿主：Rabbit

反应种属：Human

标记物：Unconjugate

克隆类型：rabbit polyclonal

技术规格

Background:	This gene belongs to the C2H2-type zinc finger gene family. The zinc finger proteins are involved in gene regulation and development, and are quite conserved throughout evolution. Like this gene product, a third of the zinc finger proteins containing C2H2 fingers also contain the KRAB domain, which has been found to be involved in protein-protein interactions.
Applications:	ELISA, IHC
Name of antibody:	ZNF43
Immunogen:	Fusion protein of human ZNF43
Full name:	zinc finger protein 43

Synonyms:	HTF6; KOX27; ZNF39L1
SwissProt:	P17038
ELISA Recommended dilution:	5000-10000
IHC positive control:	Human lung cancer and human brain
IHC Recommend dilution:	25-100



