

兔抗 PPP1R2(Ab-120/121) 多克隆抗体

中文名称：兔抗 PPP1R2(Ab-120/121) 多克隆抗体

英文名称：Anti-PPP1R2(Ab-120/121) rabbit polyclonal antibody

别名：IPP2; IPP-2

相关类别：一抗

储存：冷冻（-20℃）避光

宿主：Rabbit

抗原：PPP1R2(Ab-120/121)

反应种属：Human Mouse Rat

标记物：Unconjugate

克隆类型：Unconjugate

技术规格

Background:	Inhibitor of protein-phosphatase 1.
Applications:	WB, IHC
Name of antibody:	PPP1R2(Ab-120/121)
Immunogen:	Synthesized non-phosphopeptide derived from human PP P1R2 around the phosphorylation site of serine 120 and serine 121 (Q-E-S(p)-S(p)-G-E).
Full name:	protein phosphatase 1, regulatory (inhibitor) subunit 2
Synonyms :	IPP2; IPP-2
SwissProt:	P41236

IHC positive control:	Human breast carcinoma tissue
IHC Recommend dilution:	50-100
WB Predicted band size:	23 kDa
WB Positive control:	Jurkat cells lysate
WB Recommended dilution:	500-3000



