

兔抗 PPP2R1A 多克隆抗体

中文名称：兔抗 PPP2R1A 多克隆抗体

英文名称：Anti-PPP2R1A rabbit polyclonal antibody

别 名：protein phosphatase 2 scaffold subunit Aalpha; MRD36; PP2AA; PR65A; PP2AAALPHA; PP2A-Aalpha

相关类别：一抗

储 存：冷冻（-20℃）

宿 主：Rabbit

抗 原：PPP2R1A

反应种属：Human, Mouse

标记物：Unconjugate

克隆类型：rabbit polyclonal

技术规格

Background:	This gene encodes a constant regulatory subunit of protein phosphatase 2. Protein phosphatase 2 is one of the four major Ser/Thr phosphatases, and it is implicated in the negative control of cell growth and division. It consists of a common heteromeric core enzyme, which is composed of a catalytic subunit and a constant regulatory subunit, that associates with a variety of regulatory subunits. The constant regulatory subunit A serves as a scaffolding molecule to coordinate the assembly of the catalytic subunit.
--------------------	--

	ytic subunit and a variable regulatory B subunit. This gene encodes an alpha isoform of the constant regulatory subunit A. Alternatively spliced transcript variants have been described.
Applications:	ELISA, WB, IHC
Name of antibody:	PPP2R1A
Immunogen:	Fusion protein of human PPP2R1A
Full name:	protein phosphatase 2 scaffold subunit Aalpha
Synonyms:	MRD36; PP2AA; PR65A; PP2AAALPHA; PP2A-Aalpha
SwissProt:	P30153
ELISA Recommended dilution:	5000-10000
IHC positive control:	Human liver cancer
IHC Recommend dilution:	70-350
WB Predicted band size:	65 kDa
WB Positive control:	NIH/3T3 cell lysate
WB Recommended dilution:	1000-5000



