

兔抗 PRDX6 多克隆抗体

中文名称：兔抗 PRDX6 多克隆抗体

英文名称：Anti-PRDX6 rabbit polyclonal antibody

别名：PRX, p29, AOP2, 1-Cys, NSGPx, aiPLA2

抗原：PRDX6

储存：冷冻（-20℃）

宿主：Rabbit

相关类别：一抗

反应种属：Human, Mouse, Rat

标记物：Unconjugate

克隆类型：rabbit polyclonal

技术规格

Background:	The protein encoded by this gene is a member of the thiol-specific antioxidant protein family. This protein is a bifunctional enzyme with two distinct active sites. It is involved in redox regulation of the cell; it can reduce H ₂ O ₂ and short chain organic, fatty acid, and phospholipid hydroperoxides. It may play a role in the regulation of phospholipid turnover as well as in protection against oxidative injury.
Applications:	ELISA, WB, IHC
Name of antibody:	PRDX6
Immunogen:	Fusion protein of human PRDX6

Full name:	peroxiredoxin 6
Synonyms :	PRX, p29, AOP2, 1-Cys, NSGPx, aiPLA2
SwissProt:	P30041
ELISA Recommended dilution:	1000-5000
IHC positive control:	Human tonsil and human liver cancer
IHC Recommend dilution:	25-100
WB Predicted band size:	25 kDa
WB Positive control:	Hela, 293T, 231 and A549 cells
WB Recommended dilution:	500-2000



