

兔抗 PRKACAB 多克隆抗体

中文名称：兔抗 PRKACAB 多克隆抗体

英文名称：Anti-PRKACA/B rabbit polyclonal antibody

别名：PKACA/PKACB; PKA C-beta

相关类别：一抗

抗原：PRKACA/B

储存：冷冻（-20℃）避光

宿主：Rabbit

反应种属：Human Mouse Rat

标记物：Unconjugate

克隆类型：Unconjugate

技术规格

Background:	Phosphorylates a large number of substrates in the cytoplasm and the nucleus. Regulates the abundance of compartmentalized pools of its regulatory subunits through phosphorylation of PJA2 which binds and ubiquitinates these subunits, leading to their subsequent proteolysis. Phosphorylates CDC25B, ABL1, NFkB1, CLDN3, PSMC5/RPT6, PJA2, RYR2, RORA and VASP. RORA is activated by phosphorylation. Required for glucose-mediated adipogenic differentiation increase and osteogenic differentiation inhibition from osteoblasts. Involved in the regulation of platelets in response to thrombin and collagen; maintains circulating platelets in a resting
--------------------	---

	<p>state by phosphorylating proteins in numerous platelet inhibitory pathways when in complex with NF-kappa-B (NFKB1 and NFKB2) and I-kappa-B-alpha (NFKBIA), but thrombin and collagen disrupt these complexes and free active PRKACA stimulates platelets and leads to platelet aggregation by phosphorylating VASP. Prevents the antiproliferative and anti-invasive effects of alpha-difluoromethylornithine in breast cancer cells when activated. RYR2 channel activity is potentiated by phosphorylation in presence of luminal Ca²⁺, leading to reduced amplitude and increased frequency of store overload-induced Ca²⁺ release (SOICR) characterized by an increased rate of Ca²⁺ release and propagation velocity of spontaneous Ca²⁺ waves, despite reduced wave amplitude and resting cytosolic Ca²⁺. PSMC5/RPT6 activation by phosphorylation stimulates proteasome. Negatively regulates tight junctions (TJs) in ovarian cancer cells via CLDN3 phosphorylation. NFKB1 phosphorylation promotes NF-kappa-B p50-p50 DNA binding. Involved in embryonic development by down-regulating the Hedgehog (Hh) signaling pathway that determines embryo pattern formation and morphogenesis. Prevents meiosis resumption in prophase-arrested oocytes via CDC25B inactivation by phosphorylation. May also regulate rapid eye movement (REM) sleep in the pedunculopontine tegmental (PPT). Phosphorylates APOBEC3G and AICDA. Isoform 2 phosphorylates and activates ABL1 in sperm flagellum to promote spermatozoa capacitation.</p>
Applications:	WB
Name of antibody:	PRKACA/B
Immunogen:	Synthesized peptide derived from N-terminal of human KAPC A/B.
Full name:	protein kinase, cAMP-dependent, catalytic, alpha/beta
Synonyms :	PKACA/PKACB; PKA C-beta
SwissProt:	P17612/P22694
WB Predicted band size:	41 kDa
WB Positive control:	Jurkat cells and COLO cells lysates
WB Recommended dilution:	500-3000

