

兔抗 PRKY 多克隆抗体

中文名称：兔抗 PRKY 多克隆抗体

英文名称：Anti-PRKY rabbit polyclonal antibody

别名：PRKYP; PRKXP3

相关类别：一抗

储存：冷冻（-20℃）避光

宿主：Rabbit

抗原：PRKY

反应种属：Human

标记物：Unconjugate

克隆类型：Unconjugate

技术规格

Background:	This gene is similar to the protein kinase, X-linked gene in the pseudoautosomal region of the X chromosome. The gene is classified as a transcribed pseudogene because it has lost a coding exon that results in all transcripts being candidates for nonsense-mediated decay (NMD) and unlikely to express a protein. Abnormal recombination between this gene and a related gene on chromosome X is a frequent cause of XX males and XY females.
Applications:	WB
Name of antibody:	PRKY

Immunogen:	Synthesized peptide derived from internal of human PRKY.
Full name:	protein kinase, Y-linked, pseudogene
Synonyms :	PRKYP; PRKXP3
SwissProt:	O43930
WB Predicted band size:	32 kDa
WB Positive control:	Jurkat cells lysate
WB Recommended dilution:	500-3000

