

兔抗 PRL 多克隆抗体

中文名称：兔抗 PRL 多克隆抗体

英文名称：Anti-PRL rabbit polyclonal antibody

抗 原：PRL

储 存：冷冻（-20℃）

宿 主：Rabbit

相关类别：一抗

反应种属：Human

标记物：Unconjugate

克隆类型：rabbit polyclonal

技术规格

Background:	This gene encodes the anterior pituitary hormone prolactin. This secreted hormone is a growth regulator for many tissues, including cells of the immune system. It may also play a role in cell survival by suppressing apoptosis, and it is essential for lactation. Alternative splicing results in multiple transcript variants that encode the same protein.
Applications:	ELISA, IHC
Name of antibody:	PRL
Immunogen:	Fusion protein of human PRL
Full name:	prolactin
SwissProt:	P01236
IHC positive control:	Human gastric cancer

IHC Recommend dilution:	50-200
-------------------------	--------

