

兔抗 PRMT2 多克隆抗体

中文名称：兔抗 PRMT2 多克隆抗体

英文名称：Anti-PRMT2 rabbit polyclonal antibody

别名：HRMT1L1

抗原：PRMT2

储存：冷冻（-20℃）

宿主：Rabbit

相关类别：一抗

反应种属：Human

标记物：Unconjugate

克隆类型：rabbit polyclonal

技术规格

Background:

Protein arginine N-methyltransferase 2 is an enzyme that in humans is encoded by the PRMT2 gene. Arginine methyltransferase that methylates the guanidino nitrogens of arginyl residues in proteins such as STAT3, FBL, histone H4. Acts as a coactivator (with NCOA2) of the androgen receptor (AR)-mediated transactivation. Acts as a coactivator (with estrogen) of estrogen receptor (ER)-mediated transactivation. Enhances PGR, PPARG, RARA-mediated transactivation. May inhibit NF-kappa-B transcription and promote apoptosis. Represses E2F1 transcri

	ptional activity (in a RB1-dependent manner). May be involved in growth regulation.
Applications:	ELISA, IHC
Name of antibody:	PRMT2
Immunogen:	Fusion protein of human PRMT2
Full name:	protein arginine methyltransferase 2
Synonyms :	HRMT1L1
SwissProt:	P55345
ELISA Recommended dilution:	2000-5000
IHC positive control:	Human tonsil and human ovarian cancer
IHC Recommend dilution:	25-100

