

## 兔抗 PRPF39 多克隆抗体

中文名称: 兔抗 PRPF39 多克隆抗体

英文名称: Anti-PRPF39 rabbit polyclonal antibody

储 存: 冷冻 (-20℃) 避光

抗 原: PRPF39

宿 主: Rabbit

反应种属: Human

相关类别: 一抗

标 记 物: Unconjugate

克隆类型: rabbit polyclonal

### 技术规格

<b>Background:</b>	Involved in pre-mRNA splicing By similarity.
<b>Applications:</b>	WB
<b>Name of antibody:</b>	PRPF39
<b>Immunogen:</b>	Synthesized peptide derived from internal of human PRPF39.
<b>Full name:</b>	pre-mRNA processing factor 39
<b>Synonyms :</b>	0
<b>SwissProt:</b>	Q86UA1
<b>WB Predicted band size:</b>	78 kDa
<b>WB Positive control:</b>	K562 cells lysate
<b>WB Recommended dilution:</b>	500-3000

