

兔抗 PRPSAP1 多克隆抗体

中文名称：兔抗 PRPSAP1 多克隆抗体

英文名称：Anti-PRPSAP1 rabbit polyclonal antibody

别名：PAP39

相关类别：一抗

储存：冷冻（-20℃）

抗原：PRPSAP1

宿主：Rabbit

反应种属：Human, Mouse, Rat

标记物：Unconjugate

克隆类型：rabbit polyclonal

技术规格

Background:	PRPSAP1 (phosphoribosyl pyrophosphate synthetase-associated protein 1), also known as PAP39, is a 356 amino acid ubiquitous protein belonging to the ribose-5-phosphate pyrophosphokinase family. PRPSAP1 may play a negative regulatory role in 5-phosphoribose 1-diphosphate synthesis and is encoded by a gene mapping to human chromosome 17q25.1.
Applications:	ELISA, IHC
Name of antibody:	PRPSAP1

Immunogen:	Fusion protein of human PRPSAP1
Full name:	phosphoribosyl pyrophosphate synthetase-associated protein 1
Synonyms :	PAP39
SwissProt:	Q14558
ELISA Recommended dilution:	2000-10000
IHC positive control:	Human thyroid cancer and human liver cancer
IHC Recommend dilution:	100-300

