

兔抗 PTGIS 多克隆抗体

中文名称: 兔抗 PTGIS 多克隆抗体

英文名称: Anti-PTGIS rabbit polyclonal antibody

别 名: CYP8; PGIS; PTGI; CYP8A1

相关类别: 一抗

储 存: 冷冻(-20℃)

宿 主: Rabbit

抗 原: PTGIS

反应种属: Human, Mouse

标记物: Unconjugate

克隆类型: rabbit polyclonal

技术规格

Background:	This gene encodes a member of the cytochrome P450 s uperfamily of enzymes. The cytochrome P450 proteins ar e monooxygenases which catalyze many reactions involved in drug metabolism and synthesis of cholesterol, steroids and other lipids. However, this protein is considered a member of the cytochrome P450 superfamily on the b asis of sequence similarity rather than functional similarity.
Applications:	ELISA, WB, IHC
Name of antibody:	PTGIS



Immunogen:	Synthetic peptide of human PTGIS
Full name:	prostaglandin I2 (prostacyclin) synthase
Synonyms:	CYP8; PGIS; PTGI; CYP8A1
SwissProt:	Q16647
ELISA Recommended dilution:	5000-10000
IHC positive control:	Human cervical cancer and human ovarian cancer
IHC Recommend dilution:	20-100
WB Predicted band size:	57 kDa
WB Positive control:	Mouse lung tissue
WB Recommended dilution:	200-1000







