

兔抗 PTH1R 多克隆抗体

中文名称：兔抗 PTH1R 多克隆抗体

英文名称：Anti-PTH1R rabbit polyclonal antibody

别名：parathyroid hormone 1 receptor; PFE; PTHR; PTHR1

相关类别：一抗

储存：冷冻（-20℃）

宿主：Rabbit

抗原：PTH1R

反应种属：Human, Mouse, Rat

标记物：Unconjugate

克隆类型：rabbit polyclonal

技术规格

Background:

The protein encoded by this gene is a member of the G-protein coupled receptor family 2. This protein is a receptor for parathyroid hormone (PTH) and for parathyroid hormone-like hormone (PTHrP). The activity of this receptor is mediated by G proteins which activate adenylyl cyclase and also a phosphatidylinositol-calcium second messenger system. Defects in this receptor are known to be the cause of Jansen's metaphyseal chondrodysplasia (JMC), chondrodysplasia Blomstrand type (BOCD), as well as enchondromatosis. Two transcript variants encoding the same protein

	have been found for this gene.
Applications:	ELISA, IHC
Name of antibody:	PTH1R
Immunogen:	Fusion protein of human PTH1R
Full name:	parathyroid hormone 1 receptor
Synonyms:	PFE; PTHR; PTHR1
SwissProt:	Q03431
ELISA Recommended dilution:	5000-10000
IHC positive control:	Human tonsil and Human esophagus cancer
IHC Recommend dilution:	50-300



