

兔抗 L1RE1 多克隆抗体

中文名称：兔抗 L1RE1 多克隆抗体

英文名称：Anti-L1RE1 rabbit polyclonal antibody

别名：LINE1 retrotransposable element 1; L1.2; LRE1

相关类别：一抗

储存：冷冻（-20℃）

宿主：Rabbit

抗原：L1RE1

反应种属：Human

标记物：Unconjugate

克隆类型：rabbit polyclonal

技术规格

Background:	Nucleic acid-binding protein which is essential for retrotransposition of LINE-1 elements in the genome. May function as a nucleic acid chaperone binding its own transcript and therefore preferentially mobilizing the transcript from which they are encoded.
Applications:	ELISA, WB, IHC
Name of antibody:	L1RE1
Immunogen:	Synthetic peptide of human L1RE1
Full name:	LINE1 retrotransposable element 1
Synonyms:	L1.2; LRE1

SwissProt:	Q9UN81
ELISA Recommended dilution:	5000-10000
IHC positive control:	Human liver cancer and human ovarian cancer
IHC Recommend dilution :	25-100
WB Predicted band size:	40 kDa
WB Positive control:	A431 and 293T cell lysates
WB Recommended dilution:	200-1000





