

## 兔抗 LAIR2 多克隆抗体

中文名称：兔抗 LAIR2 多克隆抗体

英文名称：Anti-LAIR2 rabbit polyclonal antibody

别名：CD306

相关类别：一抗

储存：冷冻（-20℃）

宿主：Rabbit

抗原：LAIR2

反应种属：Human

标记物：Unconjugate

克隆类型：rabbit polyclonal

### 技术规格

<b>Background:</b>	The protein encoded by this gene is a member of the immunoglobulin superfamily. It was identified by its similarity to LAIR1, an inhibitory receptor present on mononuclear leukocytes. This gene maps to a region of 19q13.4, termed the leukocyte receptor cluster, which contains 29 genes in the immunoglobulin superfamily, including LAIR1. The function of this protein is unknown, although it is thought to be secreted and may help modulate mucosal tolerance. Two transcript variants encoding different isoforms have been found for this gene.
<b>Applications:</b>	ELISA, IHC
<b>Name of antibody:</b>	LAIR2

<b>Immunogen:</b>	Synthetic peptide of human LAIR2
<b>Full name:</b>	leukocyte-associated immunoglobulin-like receptor 2
<b>Synonyms :</b>	CD306
<b>SwissProt:</b>	Q6ISS4
<b>ELISA Recommended dilution:</b>	1000-2000
<b>IHC positive control:</b>	Human tonsil and Human brain
<b>IHC Recommend dilution:</b>	25-100

