

兔抗 ZNF581 多克隆抗体

中文名称：兔抗 ZNF581 多克隆抗体

英文名称：Anti-ZNF581 rabbit polyclonal antibody

别名：zinc finger protein 581; HSPC189

相关类别：一抗

储存：冷冻（-20℃）

宿主：Rabbit

抗原：ZNF581

反应种属：Human

标记物：Unconjugate

克隆类型：rabbit polyclonal

技术规格

Background:	ZNF581 contains 4 C2H2-type zinc fingers. ZNF581 may be involved in transcriptional regulation.
Applications:	ELISA, IHC
Name of antibody:	ZNF581
Immunogen:	Fusion protein of human ZNF581
Full name:	zinc finger protein 581
Synonyms:	HSPC189
SwissProt:	Q9P0T4
ELISA Recommended dilution:	5000-10000
IHC positive control:	Human tonsil

IHC Recommend dilution:

50-300

