

兔抗 ZNF682 多克隆抗体

中文名称：兔抗 ZNF682 多克隆抗体

英文名称：Anti-ZNF682 rabbit polyclonal antibody

别 名：BC39498_3

相关类别：一抗

储 存：冷冻（-20℃）避光

抗 原：ZNF682

宿 主：Rabbit

反应种属：Human

标记物：Unconjugate

克隆类型：Unconjugate

技术规格

Background:	May be involved in transcriptional regulation.
Applications:	WB
Name of antibody:	ZNF682
Immunogen:	Synthesized peptide derived from internal of human ZNF682 .
Full name:	zinc finger protein 682
Synonyms :	BC39498_3
SwissProt:	O95780
WB Predicted band size:	58 kDa
WB Positive control:	COLO205 and Huvec cells lysates
WB Recommended dilution:	500-3000

