

兔抗 TRIM59 多克隆抗体

中文名称：兔抗 TRIM59 多克隆抗体

英文名称：Anti-TRIM59 rabbit polyclonal antibody

别名：MRF1; TSBF1; RNF104; TRIM57

相关类别：一抗

储存：冷冻（-20℃）避光

抗原：TRIM59

宿主：Rabbit

反应种属：Human Mouse

标记物：Unconjugate

克隆类型：Unconjugate

技术规格

Background:	Belongs to the TRIM/RBCC family. Note: This description may include information from UniProtKB.
Applications:	WB, IHC, IF
Name of antibody:	TRIM59
Immunogen:	Synthesized peptide derived from internal of human TRIM59.
Full name:	tripartite motif containing 59
Synonyms :	MRF1; TSBF1; RNF104; TRIM57
SwissProt:	Q8IWR1

IHC positive control:	Human brain tissue
IHC Recommend dilution:	50-100
WB Predicted band size:	47 kDa
WB Positive control:	NIH/3T3 cells lysate
WB Recommended dilution:	500-3000
IF positive control:	Hela cells
IF Recommend dilution:	100-500



