

兔抗 LBR 多克隆抗体

中文名称: 兔抗 LBR 多克隆抗体

英文名称: Anti-LBR rabbit polyclonal antibody

别 名: PHA; LMN2R; TDRD18; DHCR14B

相关类别: 一抗

储 存: 冷冻(-20℃)

宿 主: Rabbit

抗 原: LBR

反应种属: Human

标 记 物: Unconjugate

克隆类型: rabbit polyclonal

技术规格

Background:	The protein encoded by this gene belongs to the ERG4/ER G24 family. It localized in the nuclear envelope inner mem brane and anchors the lamina and the heterochromatin to the membrane. It may mediate interaction between chroma tin and lamin B. Mutations of this gene has been associate d with autosomal recessive HEM/Greenberg skeletal dysplas ia. Alternative splicing occurs at this locus and two transcript variants encoding the same protein have been identified .
Applications:	elisa, wb, ihc



Name of antibody:	LBR
Immunogen:	Fusion protein of human LBR
Full name:	lamin B receptor
Synonyms:	PHA; LMN2R; TDRD18; DHCR14B
SwissProt:	Q14739
ELISA Recommended dilutio n:	2000-5000
IHC positive control:	Human ovarian cancer
IHC Recommend dilution:	50-200
WB Predicted band size:	71 kDa
WB Positive control:	Human fetal brain tissue
WB Recommended dilution:	500-2000





