

兔抗 ZSCAN16 多克隆抗体

中文名称：兔抗 ZSCAN16 多克隆抗体

英文名称：Anti-ZSCAN16 rabbit polyclonal antibody

别名：ZNF392; ZNF435; dJ265C24.3

抗原：ZSCAN16

储存：冷冻（-20℃）避光

宿主：Rabbit

相关类别：一抗

反应种属：Human Mouse

标记物：Unconjugate

克隆类型：rabbit polyclonal

技术规格

Background:	May be involved in transcriptional regulation.
Applications:	WB
Name of antibody:	ZSCAN16
Immunogen:	Synthesized peptide derived from internal of human ZNF435.
Full name:	zinc finger and SCAN domain containing 16
Synonyms :	ZNF392; ZNF435; dJ265C24.3
SwissProt:	Q9H4T2
WB Predicted band size:	41 kDa
WB Positive control:	NIH/3T3 cells and HT29 cells lysates
WB Recommended dilution:	500-3000

