

兔抗 ZSCAN22 多克隆抗体

中文名称：兔抗 ZSCAN22 多克隆抗体

英文名称：Anti-ZSCAN22 rabbit polyclonal antibody

相关类别：一抗

储 存：冷冻（-20℃）避光

宿 主：Rabbit

抗 原：ZSCAN22

反应种属：Human

标 记 物：Unconjugate

克隆类型：Unconjugate

技术规格

Background:	May be involved in transcriptional regulation.
Applications:	WB, IHC
Name of antibody:	ZSCAN22
Immunogen:	Synthesized peptide derived from internal of human ZSCAN22.
Full name:	hypothetical protein
Synonyms :	0
SwissProt:	P10073
IHC positive control:	Human breast carcinoma tissue
IHC Recommend dilution:	50-100
WB Predicted band size:	55 kDa
WB Positive control:	K562 cells lysate
WB Recommended dilution:	500-3000



