

兔抗 LGALS2 多克隆抗体

中文名称：	兔抗 LGALS2 多克隆抗体
英文名称：	Anti-LGALS2 rabbit polyclonal antibody
别 名：	HL14
抗 原：	LGALS2
储 存：	冷冻（-20℃）
宿 主：	Rabbit
相关类别：	一抗
反应种属：	Human
标记物：	Unconjugate
克隆类型：	rabbit polyclonal

技术规格

Background:	The protein encoded by this gene is a soluble beta-galactoside binding lectin. The encoded protein is found as a homodimer and can bind to lymphotoxin-alpha. A single nucleotide polymorphism in an intron of this gene can alter the transcriptional level of the protein, with a resultant increased risk of myocardial infarction.
Applications:	ELISA, WB, IHC
Name of antibody:	LGALS2
Immunogen:	Fusion protein of human LGALS2
Full name:	galectin 2
Synonyms:	HL14

SwissProt:	P05162
ELISA Recommended dilution:	5000-10000
IHC positive control:	Human liver cancer
IHC Recommend dilution:	50-200
WB Predicted band size:	15 kDa
WB Positive control:	Human fetal intestines tissue lysate
WB Recommended dilution:	500-2000



