

## 小鼠抗 ACTR3 单克隆抗体

中文名称：小鼠抗 ACTR3 单克隆抗体

英文名称：Anti-ACTR3 mouse monoclonal antibody

别名：ARP3 actin related protein 3 homolog; ARP3

相关类别：一抗

储存：冷冻（-20℃）

宿主：Mouse

抗原：ACTR3

反应种属：Human

标记物：Unconjugate

克隆类型：mouse monoclonal

### 技术规格

<b>Background:</b>	The specific function of this gene has not yet been determined; however, the protein it encodes is known to be a major constituent of the ARP2/3 complex. This complex is located at the cell surface and is essential to cell shape and motility through lamellipodial actin assembly and protrusion. Three transcript variants encoding two different isoforms have been found for this gene.
<b>Applications:</b>	WB, IHC
<b>Name of antibody:</b>	ACTR3
<b>Immunogen:</b>	Fusion protein of human ACTR3
<b>Full name:</b>	ARP3 actin related protein 3 homolog
<b>Synonyms:</b>	ARP3
<b>SwissProt:</b>	P61158

<b>IHC positive control:</b>	Human tonsil tissue and Human spleen tissue
<b>IHC Recommend dilution:</b>	500-2000
<b>WB Predicted band size:</b>	47 KD
<b>WB Positive control:</b>	Hela, K562, A431 cell lysates
<b>WB Recommended dilution:</b>	500-2000