

## 小鼠抗 CHAF1B 单克隆抗体

中文名称: 小鼠抗 CHAF1B 单克隆抗体

英文名称: Anti-CHAF1B mouse monoclonal antibody

别名: CAF-1; CAF-IP60; CAF1; CAF1A; CAF1P60; MPHOSPH7; MPP7

抗原: CHAF1B

储存: 冷冻 (-20℃) 避光

宿主: Mouse

反应种属: Human , Dog , Rat , Monkey , Mouse

相关类别: 一抗

标记物: Unconjugate

### 技术规格

<b>Background:</b>	Chromatin assembly factor I (CAF-I) is required for the assembly of histone octamers onto newly-replicated DNA. CAF-I is composed of three protein subunits, p50, p60, and p150. The protein encoded by this gene corresponds to the p60 subunit and is required for chromatin assembly after replication. The encoded protein is differentially phosphorylated in a cell cycle-dependent manner. In addition, it is normally found in the nucleus except during mitosis, when it is released into the cytoplasm. This protein is a member of the WD-repeat HIR1 family and may also be involved in DNA repair. [provided by RefSeq, Jul 2008].
<b>Applications:</b>	WB, IHC, IF

<b>Name of antibody:</b>	CHAF1B
<b>Immunogen:</b>	Fusion protein of human CHAF1B
<b>Full name:</b>	chromatin assembly factor 1, subunit B (p60) (CHAF1B)
<b>Synonyms:</b>	CAF-1; CAF-IP60; CAF1; CAF1A; CAF1P60; MPHOSPH7; MP P7
<b>SwissProt:</b>	Q13112
<b>IHC positive control:</b>	carcinoma of human prostate tissue and human lymph node tissue; human lymphoma tissue
<b>IHC Recommend dilution:</b>	30-150
<b>WB Predicted band size:</b>	61 kDa
<b>WB Positive control:</b>	HepG2, Hela, SVT2, A549, COS7, Jurkat, MDCK, PC12, MC F-7 cell lysates
<b>WB Recommended dilution:</b>	200-1000