

## 小鼠抗 CITED1 单克隆抗体

中文名称： 小鼠抗 CITED1 单克隆抗体

英文名称： Anti-CITED1 mouse monoclonal antibody

别 名： MSG1

相关类别： 一抗

储 存： 冷冻（-20℃） 避光

宿 主： Mouse

抗 原： CITED1

反应种属： Human

标 记 物： Unconjugate

克隆类型： mouse monoclonal

### 技术规格

<b>Background:</b>	This gene encodes a member of the CREB-binding protein/p300-interacting transactivator with Asp/Glu-rich C-terminal domain (CITED) family of proteins. The encoded protein, also known as melanocyte-specific gene 1, may function as a transcriptional coactivator and may play a role in pigmentation of melanocytes. Alternatively spliced transcript variants have been described. [provided by RefSeq, Jan 2009].
<b>Applications:</b>	WB, IHC
<b>Name of antibody:</b>	CITED1
<b>Immunogen:</b>	Fusion protein of human CITED1

<b>Full name:</b>	Cbp/p300-interacting transactivator, with Glu/Asp-rich carboxy-terminal domain, 1 (CITED1), transcript variant 1
<b>Synonyms:</b>	MSG1
<b>SwissProt:</b>	Q99966
<b>IHC positive control:</b>	human testicular cancer tissue
<b>IHC Recommend dilution:</b>	30-150
<b>WB Predicted band size:</b>	20 kDa
<b>WB Positive control:</b>	Hela cell lysate
<b>WB Recommended dilution:</b>	200-1000