

小鼠抗 DNAJB4 单克隆抗体

中文名称:小鼠抗 DNAJB4 单克隆抗体

英文名称: Anti-DNAJB4 mouse monoclonal antibody

别 名: DnaJ heat shock protein family (Hsp40) member B4; DjB4; HLJ1; DNAJW

相关类别: 一抗

储 存: 冷冻(-20℃)

宿 主: Mouse

抗原: DNAJB4

反应种属: Human

标记物: Unconjugate

克隆类型: mouse monoclonal

技术规格

| Background: | The protein encoded by this gene is a molecular chapero ne, tumor suppressor, and member of the heat shock pro tein-40 family. The encoded protein binds the cell adhesi on protein E-cadherin and targets it to the plasma memb rane. This protein also binds incorrectly folded E-cadherin and targets it for endoplasmic reticulum-associated degra dation. This gene is a strong tumor suppressor for colore ctal carcinoma, and downregulation of it may serve as a good biomarker for predicting patient outcomes. Several transcript variants encoding different isoforms have been found for this gene. |
|-------------------|---|
| Applications: | WB |
| Name of antibody: | DNAJB4 |



| Immunogen: | Fusion protein of human DNAJB4 |
|--------------------------|--|
| Full name: | DnaJ heat shock protein family (Hsp40) member B4 |
| Synonyms: | DjB4; HLJ1; DNAJW |
| SwissProt: | Q9UDY4 |
| WB Predicted band size: | 38 KD |
| WB Positive control: | HepG2 cell lysate |
| WB Recommended dilution: | 500-2000 |