

小鼠抗 ENO1 单克隆抗体

中文名称： 小鼠抗 ENO1 单克隆抗体

英文名称： Anti-ENO1 mouse monoclonal antibody

别名： enolase 1; NNE; PPH; MPB1; ENO1L1; HEL-S-17

相关类别： 一抗

储存： 冷冻（-20℃）

宿主： Mouse

抗原： ENO1

反应种属： Human, Mouse, Rat, Monkey

标记物： Unconjugate

克隆类型： Mouse Monoclonal

技术规格

Background:	This gene encodes alpha-enolase, one of three enolase isoenzymes found in mammals. Each isoenzyme is a homodimer composed of 2 alpha, 2 gamma, or 2 beta subunits, and functions as a glycolytic enzyme. Alpha-enolase in addition, functions as a structural lens protein (tau-crystallin) in the monomeric form. Alternative splicing of this gene results in a shorter isoform that has been shown to bind to the c-myc promoter and function as a tumor suppressor. Several pseudogenes have been identified, including one on the long arm of chromosome 1. Alpha-enolase has also been identified as an autoantigen i
--------------------	---

	n Hashimoto encephalopathy.
Applications:	WB
Name of antibody:	ENO1
Immunogen:	Fusion protein of human ENO1
Full name:	enolase 1
Synonyms:	NNE; PPH; MPB1; ENO1L1; HEL-S-17
SwissProt:	P06733
WB Predicted band size:	47 kDa
WB Positive control:	T47D, 3T3, COS7, C6 and Hela cell lysates
WB Recommended dilution:	1000-5000

