

兔抗 LSP1 多克隆抗体

中文名称：兔抗 LSP1 多克隆抗体

英文名称：Anti-LSP1 rabbit polyclonal antibody

别名：WP34; pp52

抗原：LSP1

储存：冷冻（-20℃）避光

宿主：Rabbit

反应种属：Human

相关类别：一抗

标记物：Unconjugate

克隆类型：Unconjugate

技术规格

Background:	This gene encodes an intracellular F-actin binding protein. The protein is expressed in lymphocytes, neutrophils, macrophages, and endothelium and may regulate neutrophil motility, adhesion to fibrinogen matrix proteins, and transendothelial migration. Alternative splicing results in multiple transcript variants encoding different isoforms.
Applications:	WB
Name of antibody:	LSP1
Immunogen:	Synthesized peptide derived from internal of human LSP1.
Full name:	lymphocyte-specific protein 1

Synonyms :	WP34; pp52
SwissProt:	P33241
WB Predicted band size:	37 kDa
WB Positive control:	Huvec cells and K562 cells lysates
WB Recommended dilution:	500-3000

