

## 小鼠抗 PDCL3 单克隆抗体

中文名称: 小鼠抗 PDCL3 单克隆抗体

英文名称: Anti-PDCL3 mouse monoclonal antibody

别 名: phosducin like 3; VIAF; PHLP3; VIAF1; HTPHLP; PHLP2A

相关类别: 一抗

储 存: 冷冻(-20℃)

宿 主: Mouse

抗 原: PDCL3

反应种属: Human

标 记 物: Unconjugate

克隆类型: mouse monoclonal

## 技术规格

Background:	This gene encodes a member of the phosducin-like protein f amily and is a putative modulator of heterotrimeric G proteins. The protein shares extensive amino acid sequence homology with phosducin. Members of the phosducin-like protein family have been shown to bind to the beta-gamma subunit s of G proteins.
Applications:	WB
Name of antibody:	PDCL3
Immunogen:	Fusion protein of human PDCL3
Full name:	phosducin like 3
Synonyms:	VIAF; PHLP3; VIAF1; HTPHLP; PHLP2A



SwissProt:	Q9H2J4
WB Predicted band size:	28 KD
WB Positive control:	Jurkat, A549, K562, 293T, Hela cell lysates
WB Recommended dilution :	500-2000